

Dry Block Heaters 2 Block Dry Block Heaters with Lid







Multi-purpose Dry Block Heaters are ideal for applications that require temperature stability. The close tube-to-block contact enables maximum heat retention, resulting in efficient heating. High-wattage, constant-temp. analog models are an economical option, while digital models offer exceptional temp. uniformity & stability for applications that require repeatable results. Units hold optional interchangeable modular blocks with over 40 options.

- Delivering Exceptional Temperature Stability and Uniformity
- Accurate and Fast Sample Heating with PID Microprocessor Temperature Control (Digital models)
- Digital Models can be Calibrated to an External Temperature Device

- Exceptional temperature uniformity and stability
- Optional external temperature probe
- Holds interchangeable modular blocks

Designed for applications that require repeatable results and superior temperature stability. These multi-purpose units are ideal for incubation and activation of cultures, enzyme reactions, immunoassays, melting/boiling points, and a wide variety of other laboratory procedures. Integral support rod holder with locking knob accepts optional External Temperature Probe Kit. Optional External Temperature Probe Kit monitors actual block or sample temperature. Each of the five models accepts separate interchangeable modular blocks, accommodating various tube sizes from 0.2 mL micro-tubes to 50 mL centrifuge tubes. Each block has a thermometer well for measuring block temperature. Efficient heating due to close tube-and-block contact for maximum heat retention. Heaters require OHAUS modular heating blocks for operation. Modular blocks are sold separately (see pages 73–75).

### **Operating Features:**

**Microprocessor Control:** PID temperature control, with optional external RTD probe, offers a temperature stability as low as  $\pm$  0.1°C with a temperature uniformity as low as  $\pm$  0.1°C. Samples are heated to temperature quickly and accurately. Temperature is adjusted in  $\pm$  0.1°C increments.

**LED Display:** Touch pad controls with easy-to-read, independent LED displays for temperature, and time. Provides repeatable and accurate results every time and is easily visible across lab benches. Shows set-point and actual temperature. Timer will display elapsed time or, when programmed to user defined limit, will shut off unit when time reaches zero. Display will show last used settings, even after power has been turned off.

**Temperature Calibration Mode:** Allows user to calibrate unit to an external temperature device.

#### **Safety Features:**

**Caution Hot Indicator:** Hot warning symbol light is illuminated when the temperature is above 40°C.

**Overshoot Protection:** If the unit exceeds the set temperature by 10°C the unit will automatically stop heating.

**Audible Alarm:** In timed mode, alarm will sound when time reaches zero and when unit reaches set-point temperature.

**CAUTION!** To avoid possible electrical hazard, do not fill well or block opening with water or other fluids. Unit is designed as a dry bath/incubator.

### **Operating Conditions:**

Units can be run in conditions from 18 to 33°C (64 to 91°F), 20% to 80% relative humidity, non-condensing.



### **Applications:**

Denaturing proteins, DNA applications, ELISA and other immunoassay studies.

#### **Ordering Information:**

Units include a detachable, 3-wire cord and plug. Modular blocks are sold separately (see pages 73-75).

Size	Temperature Range	Temperature Stability @ 37°C	Uniformity Within the Block @ 37°C	Uniformity Across Similar Blocks @ 37°C	Heat-Up Time to 100°C
1 block				N/A	45 minutes
2 block	Ambient +5°C	Ambient $+5^{\circ}$ C 120 volt units: to 120°C $\pm 0.1^{\circ}$ C	120 volt units:	± 0.1°C	50 minutes
4 block	to 120°C		± 0.1°C	± 0.2°C	60 minutes
6 block				± 0.3°C	65 minutes

## **Dry Block Heaters**

### **Advanced Dry Block Heaters**

Block Capacity	Dimensions (L $\times$ W $\times$ H)	Model	Item Number	Item Number
1	12.4 × 8 × 3.5" (31.5 × 20.3 × 8.9 cm)	HB1DG	30392061	VN, MY, SG, ID, IN
1	$12.4 \times 8 \times 3.5$ " (31.5 × 20.3 × 8.9 cm)	HB1DG	30392062	AU, NZ
1	12.4 × 8 × 3.5" (31.5 × 20.3 × 8.9 cm)	HB1DG	30392064	TH, PH
2	15.4 × 8 × 3.5" (39.1 × 20.3 × 8.9 cm)	HB2DG	30392082	VN, MY, SG, ID, IN
2	$15.4 \times 8 \times 3.5$ " (39.1 × 20.3 × 8.9 cm)	HB2DG	30392083	AU, NZ
2	$15.4 \times 8 \times 3.5$ " (39.1 × 20.3 × 8.9 cm)	HB2DG	30392085	TH, PH
4	$16.9 \times 8 \times 3.5$ " (42.9 × 20.3 × 8.9 cm)	HB4DG	30392089	VN, MY, SG, ID, IN
4	$16.9 \times 8 \times 3.5$ " (42.9 × 20.3 × 8.9 cm)	HB4DG	30392090	AU, NZ
4	16.9 × 8 × 3.5" (42.9 × 20.3 × 8.9 cm)	HB4DG	30392092	TH, PH
6	$20.9 \times 8 \times 3.5$ " (53.1 × 20.3 × 8.9 cm)	HB6DG	30392096	VN, MY, SG, ID, IN
6	$20.9 \times 8 \times 3.5$ " (53.1 × 20.3 × 8.9 cm)	HB6DG	30392097	AU, NZ
6	$20.9 \times 8 \times 3.5$ " (53.1 × 20.3 × 8.9 cm)	HB6DG	30392099	TH, PH

### **External Temperature Probe Kit**

Enables the unit to read actual block or sample temperature and display that temperature on the control panel. The optional External Temperature Probe Kit includes a stainless steel RTD PT100 probe, 18" (45.7 cm) stainless steel support rod, thermometer/temperature probe extension clamp and hook connector. The PT100 RTD probe is designed to fit perfectly in to the thermometer well of each modular block.



Description	Item Number
Optional External Temperature Probe Kit	30400246



- Multi-purpose use
- Holds interchangeable modular blocks
- Analog controls

High wattage, constant-temperature Dry Block Heaters are economical, versatile and compact. These multi-purpose units are ideal for incubation and activation of cultures, enzyme reactions, immunoassays, melting/boiling points and a wide variety of other laboratory procedures. Each of the five models accept separate interchangeable modular blocks, accommodating various sample enclosures such as micro-tubes, centrifuge tubes, vials, microplates, and PCR strips or tubes. Each block has a thermometer well for measuring block temperature. Anodized aluminum modular blocks provide superior temperature stability and heat transfer. Heaters require OHAUS modular heating blocks for operation. Modular blocks are sold separately (see pages 73-75).

### **Operating Features:**

High wattage, constant-temperature Dry Block Heaters are economical, versatile and compact.

**Microprocessor Control:** PID temperature controller maintains precise temperature control. Samples are heated to temperature quickly and accurately. **Adjustment Knobs:** Dual temperature control knobs with dial markings from 1 to 10 for low temperature and high temperature adjustments. Low range knob adjusts from ambient to 100°C and high range knob adjusts from 75°C to 150°C. **CAUTION!** To avoid possible electrical hazard, do not fill well or block opening with water or other fluids. Unit is designed as a dry bath/incubator.

#### **Operating Conditions:**

Units can be run in conditions from 18 to  $33^{\circ}$ C (64 to  $91^{\circ}$ F), 20% to 80% relative humidity, non-condensing.

#### **Applications:**

Coagulation and RH Studies.

### **Ordering Information:**

Units include a detachable, 3-wire cord and plug``. Modular blocks are sold separately (see pages 73-75).



Size	Temperature Range	Uniformity Within the Block @ 37°C	Temperature Stability @ 37°C	Uniformity Across Similar Blocks @ 37°C	Heat-Up Time to 100°C	
1 block	Low Range:		120 volt units: ± 1.0°C	N/A	45 minutes	
2 block	Ambient +5°C to 100°C	120 volt units: ± 0.1°C	120 1011 011101	120 volt units: ± 1.5°C	± 0.1°C	50 minutes
4 block	High Range:			120 volt units: ± 2.0°C	± 0.2°C	70 minutes
6 block	75℃ to 150℃		120 VOIL UNITS: ± 2.0 C	± 0.3°C	75 minutes	

## **Dry Block Heaters**

### **Standard Dry Block Heaters**

Block Capacity	Dimensions (L × W × H)	Model	Item Number	Country Code
1	12.4 × 8 × 3.5" (31.5 × 20.3 × 8.9 cm)	HB1AL	30392047	VN, MY, SG, ID, IN
1	12.4 × 8 × 3.5" (31.5 × 20.3 × 8.9 cm)	HB1AL	30392048	AU, NZ
1	$12.4 \times 8 \times 3.5$ " (31.5 × 20.3 × 8.9 cm)	HB1AL	30392050	TH, PH
2	15.4 × 8 × 3.5" (39.1 × 20.3 × 8.9 cm)	HB2AL	30392054	VN, MY, SG, ID, IN
2	15.4 × 8 × 3.5" (39.1 × 20.3 × 8.9 cm)	HB2AL	30392055	AU, NZ
2	$15.4 \times 8 \times 3.5$ " (39.1 × 20.3 × 8.9 cm)	HB2AL	30392057	TH, PH
4	$16.9 \times 8 \times 3.5$ " ( $42.9 \times 20.3 \times 8.9$ cm)	HB4AL	30392068	VN, MY, SG, ID, IN
4	$16.9 \times 8 \times 3.5$ " ( $42.9 \times 20.3 \times 8.9$ cm)	HB4AL	30392069	AU, NZ
4	16.9 × 8 × 3.5" (42.9 × 20.3 × 8.9 cm)	HB4AL	30392071	TH, PH
6	20.9 × 8 × 3.5" (53.1 × 20.3 × 8.9 cm)	HB6AL	30392075	VN, MY, SG, ID, IN
6	20.9 × 8 × 3.5" (53.1 × 20.3 × 8.9 cm)	HB6AL	30392076	AU, NZ
6	20.9 × 8 × 3.5" (53.1 × 20.3 × 8.9 cm)	HB6AL	30392078	TH, PH

2 Block Dry Block Heaters with Lid



Multi-purpose Dry Block Heaters are ideal for applications that require temperature stability. The close tube-to-block contact enables maximum heat retention, resulting in efficient heating. High-wattage, constant-temp. Analog models are an economical option, while digital models offer exceptional temp. uniformity & stability for applications that require repeatable results. Units hold optional interchangeable modular blocks with over 40 options.

- Delivering Exceptional Temperature Stability and Uniformity
- Heated Lid Model Reduces Condensation on Lids to Maintain Sample Integrity
- Calibrate to an External Temperature Device with Temperature Calibration Mode

### 2 Block Dry Block Heaters with Lid

- Exceptional uniformity, stability, and regulation of temperature
- Heated lid reduces condensation on sample lids
- Optional external temperature probe

OHAUS Digital Dry Block Heaters with Heated Lid are designed for applications that require repeatable results and superior temperature stability. These multi-purpose units are ideal for isothermal incubation, enzyme reactions, immunoassays, nucleic acid denaturation and a wide variety of other laboratory procedures. Integral support rod holder with locking knob accepts optional External Temperature Probe Kit. Optional External Temperature Probe Kit monitors actual block or sample temperature. Block heater accepts one microplate block or two separate interchangeable modular blocks, accommodating various tube sizes from 0.2 mL micro-tubes to test tubes or vials up to 85 mm in height. Each block has a thermometer well for measuring block temperature. Efficient heating due to close tube-and-block contact with a low density design for exceptional temperature uniformity. The heated lid helps to regulate the temperature and reduce the amount of condensation on sample lids. Heaters require OHAUS modular heating blocks for operation. Modular blocks are sold separately (see pages 73-75).

### **Operating Features:**

**Microprocessor Control:** PID temperature control, with optional external RTD probe, offers a temperature stability as low as  $\pm$  0.1°C with a temperature uniformity as low as  $\pm$  0.1°C. Samples are heated to temperature quickly and accurately. Temperature is adjusted in  $\pm$  0.1°C increments.

**LED Display:** Touch pad controls with easy-to-read, independent LED displays for temperature, and time. Provides repeatable and accurate results every time and is easily visible across lab benches. Shows set-point and actual temperature. Timer will display elapsed time or, when programmed to user defined limit, will shut off unit when time reaches zero. Display will show last used settings, even after power has been turned off.

**Temperature Calibration Mode:** Allows user to calibrate unit to an external temperature device.

#### **Safety Features:**

**Caution Hot Indicator:** Hot warning symbol light is illuminated when the temperature is above 40°C.

**Overshoot Protection:** If the unit exceeds the set temperature by 10°C the unit will automatically stop heating.

**Audible Alarm:** In timed mode, alarm will sound when time reaches zero and when unit reaches set-point temperature.

**CAUTION!** To avoid possible electrical hazard, do not fill well or block opening with water or other fluids. Unit is designed as a dry bath/incubator.



### **Operating Conditions:**

Units can be run in environments from 18 to 33°C (64 to 91°F), 20% to 80% relative humidity, non-condensing.

### **Applications:**

Isothermal incubation, enzyme reactions, immunoassays and nucleic acid denaturation and a wide variety of other laboratory procedures.

#### **Ordering Information:**

Units include a detachable, 3-wire cord and plug. For optional temperature probe kits see page 68. Modular blocks sold separately (see pages 73–75).

Size	Temperature Range	Temperature Stability @ 37°C	Uniformity Within the Block @ 37°C	Uniformity Across Similar Blocks @ 37°C	Heat-Up Time to 100°C
2 block	Ambient +5°C to 100°C	120 volt units: ± 0.1°C	120 volt units: ± 1.0°C	± 0.1°C	50 minutes

Block Capacity	Dimensions (L × W × H)	Model	Item Number	Country Code
2	15.4 × 8 × 7" (39.1 × 20.3 × 17.8 cm)	HB2DGHL	30392103	VN, MY, SG, ID, IN
2	15.4 × 8 × 7" (39.1 × 20.3 × 17.8 cm)	HB2DGHL	30392104	AU, NZ
2	15.4 × 8 × 7" (39.1 × 20.3 × 17.8 cm)	HB2DGHL	30392106	TH, PH

### Accessories

### **Modular Blocks & Accessories**

Modular blocks are constructed from a solid anodized aluminum block\*. The close contact of tubes-to-block walls allow for maximum heat retention. Each block has a thermometer well for measuring block temperature.

**CAUTION!** To avoid possible electrical hazard, do not fill well or block with water or other fluids. Units are designed as a dry bath/incubator.

Single block dimensions (L  $\times$  W  $\times$  H): 3.75  $\times$  3  $\times$  2" (9.5  $\times$  7.6  $\times$  5.1 cm) **Double block dimensions** (L  $\times$  W  $\times$  H): 6  $\times$  3.75  $\times$  2.25" (15.2  $\times$  9.5  $\times$  5.7 cm)

### **Microcentrifuge Tube Blocks**

#### Single block

Brand/Sample Type	No. of Wells	Well Dia.	Well Depth	Item Number
0.5 mL Tube	30	7.9 mm	27.6 mm	30400157
1.5 mL Tube	20	11.1 mm	39.1 mm	30400159
1.5 mL Eppendorf™ Tube	20	11.5 mm	36.9 mm	30400162
2 mL Eppendorf™ Tube	20	11.5 mm	38.1 mm	30400191
2 mL Corning™ Tube	20	10.9 mm	38.1 mm	30400192



### **Titer Plate Block**

#### **Double block**

Fits 2/4/6 block Dry Block Heaters. Ideal for 96-well or 384-well titer plates. Recessed well for better stability, flat surface good for flat and round bottom plates.

Sample Type	Well Depth	Item Number
Titer Plate	13.5 mm	30400164




### **Conical-Bottom Centrifuge Tube Blocks**

Single/Sample Type	No. of Wells	Well Dia.	Well Depth	Item Number
15 mL Tube	12	17.1 mm	44.5 mm	30400172
50 mL Tube	5	29.0 mm	47.6 mm	30400168



### **Standard Test Tube Blocks**

Single/Sample Type	No. of Wells	Well Dia.	Well Depth	Item Number
6 mm Tube	30	8.3 mm	48.4 mm	30400158
10 mm Tube	24	10.7 mm	48.4 mm	30400151
12/13 mm Tube	20	13.9 mm	48.4 mm	30400152
12/13 mm Tube	16	13.9 mm	48.4 mm	30400165
15/16 mm Tube	12	17.5 mm	48.4 mm	30400153
17/18 mm Tube	12	19.1 mm	48.4 mm	30400195
20 mm Tube	8	21.0 mm	48.4 mm	30400154
25 mm Tube	6	26.2 mm	48.4 mm	30400155
35 mm Tube	4	35.0 mm	47.6 mm	30400167



<sup>\*</sup> Block color subject to change

### **Modular Blocks\* & Accessories**

### **Combination Blocks**

### Single block

These blocks have been designed for variable sized samples.

Sample Type		No. of Wells	Well Dia.	Well Depth	Item Number
Test Tube Combination	6 mm	6	8.3 mm	48.4 mm	
	12/13 mm	5	13.8 mm	48.4 mm	30400156
	25 mm	3	26.2 mm	48.4 mm	
Centrifuge Tube	1.5 mL	4	11.1 mm	39.1 mm	
Combination	15 mL	3	17.1 mm	44.5 mm	30400193
	50 mL	2	29.0 mm	47.6 mm	
Micro-Tube Combination	0.5 mL	6	7.9 mm	27.6 mm	
	1.5 mL	10	11.1 mm	39.1 mm	30400194
	2 mL	5	11.5 mm	38.1 mm	




### **Vial Blocks**

### Single block

Designed for sample/serum and scintillation vials.

Sample Type	No. of Wells	Well Dia.	Well Depth	Item Number
12 mm Vial	20	12.7 mm	30 mm	30400182
15 mm Vial	20	15.8 mm	35 mm	30400183
16 mm Vial	15	16.4 mm	45 mm	30400190
17 mm Vial	12	17.8 mm	45 mm	30400184
19 mm Vial	12	19.7 mm	45 mm	30400185
21 mm Vial	9	21.7 mm	45 mm	30400186
23 mm Vial	8	23.8 mm	45 mm	30400187
25 mm Vial	8	25.8 mm	45 mm	30400188
28 mm Vial	6	28.8 mm	45 mm	30400189



### PCR—Plate, Tube, Strip Blocks

### Single and double block

Tapered tube wells for 0.2 mL tubes. Spaced for easy access and removal.

Sample Type	Block Size	No. of Wells	Well Dia.	Well Depth	Item Number
96-well PCR Plate	Double	96	6.4 mm	15.5 mm	30400171
10 × 8 Tube Strips	Single	80	6.4 mm	15.5 mm	30400169
Individual Tubes	Single	64	6.4 mm	20.2 mm	30400170



### **Cuvette Block**

### Single block

Two parallel slots fit 6 cuvettes in each slot, side-by-side.

Sample Type	No. of Wells	Well Depth	Item Number
(12) 12.5 mm Cuvettes	2	25.4 mm	30400161





### **Modular Blocks\* & Accessories**

### **Solid Blocks**

### Single and double block

For use as a low-temperature hotplate, slide drying, or for custom drilling to make a custom block.

Block Size	Dimensions (L $\times$ W $\times$ H)	Item Number
Single	3.75 × 3 × 2" (9.5 × 7.6 × 5.1 cm)	30400160
Double	6 × 3.75 × 2.25" (15.2 × 9.5 × 5.7 cm)	30400166



### **Stainless Steel Sand Baths**

Ideal for irregular vessels. Stainless steel construction for superior corrosion resistance. Designed to hold sand, stainless steel shot or non-volatile fluids.

For Unit	Dimensions (L $\times$ W $\times$ H)	Item Number
1 Block Dry Block Heater	3.8 × 3.0 × 2.5" (9.5 × 7.6 × 6.4 cm)	30400173
2 Block Dry Block Heater	3.8 × 5.9 × 2.5" (9.5 × 14.9 × 6.4 cm)	30400174
4 Block Dry Block Heater	5.8 × 7.5 × 2.5" (14.7 × 19.1 × 6.4 cm)	30400175
6 Block Dry Block Heater	5.8 × 11.3 × 2.5" (14.7 × 28.7 × 6.4 cm)	30400176

Accessories	Item Number
Sand, 1 lb (0.45 kg)	30400177
Stainless Steel Shot, 1 lb (0.45 kg)	30400178



### **Low Temperature Covers**

Plexiglass cover reduces air flow for additional temperature stability in low-temperature applications. Two sides are open 0.25" (6.4 mm).

Block Size	Dimensions (L $\times$ W $\times$ H)	Item Number
2 Block Cover	6.5 × 6.5 × 1.6" (16.5 × 16.5 × 4.1 cm)	30400179
4 Block Cover	8.5 × 8.5 × 1.6" (21.6 × 21.6 × 4.1 cm)	30400180
6 Block Cover	12.5 × 8.5 × 1.6" (31.8 × 21.6 × 4.1 cm)	30400181



<sup>\*</sup> Block color subject to change

# LabJaws Clamps & Supports



Lab-Lifts

# Lab-Lifts



LabJaws offers the largest selection of Lab-Frame Kits, Lab-Lifts, Rods and Support Stands in the industry to hold glassware safely. Whether you're looking for a standard lab-frame or a customizable set-up, we have options to fit any space or application. Select from 14 preconfigured lab-frame kits (connectors and feet included) to rods which can be purchased individually and combined with LabJaws accessories to build a customized frame.

- Choose from Stainless Steel and Aluminum Options to Best Suit Your Application
- Range of Pre-Configured Lab-Frame Kits and Customizable Rods for Your Lab Space
- Lab-Lifts Provide Exceptional Stability to Hold Items at Varying Heights

## LabJaws Clamps & Supports

### Lab-Lifts

### **OHAUS Aluminum Lab-Lifts**

- · Exceptional stability and durability
- Aluminum construction
- Three convenient sizes

Aluminum Lab-Lifts provide stable height adjustment for various items in the lab such as flasks, baths, and small equipment. Top and bottom decks are constructed of anodized aluminum. Internal supports and drive screws are constructed of stainless steel. Oversized sure-grip adjustment knobs provide smooth and accurate height adjustment. Lab-Lifts accept optional Support Rod Kit which mounts to the upper deck.

Deck Size	Min. to Max. Height	Max. Load*	Item Number
4 × 4" (102 × 102 mm)	2.5 to 5" (64 to 127 mm)	66 lbs (29.94 kg)	30400007
6 × 6" (152 × 152 mm)	3 to 9.75" (76 to 248 mm)	132 lbs (59.87 kg)	30400008
8 × 8" (203 × 203 mm)	3 to 9.75" (76 to 248 mm)	176 lbs (79.83 kg)	30400009
10 × 10" (254 × 254 mm)	3.5 to 13" (89 to 330 mm)	186 lbs (84.37 kg)	30400010

<sup>\*</sup>NOTE: Maximum load rating represents static weight only. Static weight is the amount a unit can hold, not lift.









OPTIONAL ACCESSORIES

### 17" Support Rod Kit

Ideal for creating an adjustable support stand for mounting various items such as thermometer clamps, temperature probes, flask and column clamps. This kit allows you to mount a 17" (432 mm) threaded vertical support rod to the upper deck of a  $6 \times 6$ ",  $8 \times 8$ ",  $10 \times 10$ ",  $12 \times 12$ " or  $16 \times 16$ " Lab-Lift by screwing the rod into the pre-drilled hole.

### 17" Support Rod Kit includes:

- 1 17" (432 mm) Stainless Steel Rod
- 1 Jam Nut
- 1 Flat Washer

Description	Rod Diameter	Item Number
17" Support Rod Kit	0.51" (13 mm)	30400050



14

### Lab-Lifts

### **OHAUS Heavy-Duty Lab-Lifts**

- Stainless steel construction
- Seven convenient sizes to choose from
- Autoclavable and chemical resistant

These Heavy-Duty, Stainless Steel Lab-Lifts are ultra-stable lifting platforms with exceptional strength and durability. Constructed of stainless steel, Lab-Lifts are designed for use in extreme environments and high load applications. Equipped with oversized, sure-grip adjustment knobs that provide extra leverage for easy height adjustments. Durable construction allows lifts to be autoclaved or chemically cleaned. Ideal for use in fume hoods or bench tops and holds a variety of items such as glassware, hotplates, baths, and magnetic stirrers.

Deck Size	Min. to Max. Height	Max. Load*	Item Number
3 × 3" (76 × 76 mm)	2.5 to 5" (64 to 127 mm)	100 lbs (45.36 kg)	30400000
4 × 4" (102 × 102 mm)	2.5 to 5" (64 to 127 mm)	100 lbs (45.36 kg)	30400001
6 × 6" (152 × 152 mm)	3 to 9.75" (76 to 248 mm)	133 lbs (60.33 kg)	30400002
8 × 8" (203 × 203 mm)	3 to 9.75" (76 to 248 mm)	227 lbs (102.97 kg)	30400003
10 × 10" (254 × 254 mm)	3.5 to 13" (89 to 330 mm)	247 lbs (112.04 kg)	30400004
12 × 12" (305 × 305 mm)	4 to 19.5" (102 to 495 mm)	100 lbs (45.36 kg)	30400005
16 × 16" (406 × 406 mm)	4 to 19.5" (102 to 495 mm)	100 lbs (45.36 kg)	30400006

<sup>\*</sup>NOTE: Maximum load rating represents static weight only. Static weight is the amount a unit can hold, not lift.















#### OPTIONAL ACCESSORIES

#### 17" Support Rod Kit

Ideal for creating an adjustable support stand for mounting various items such as thermometer clamps, temperature probes, flask and column clamps. This kit allows you to mount a 17" (432 mm) threaded vertical support rod to the upper deck of a  $6\times6$ ",  $8\times8$ ",  $10\times10$ ",  $12\times12$ " or  $16\times16$ " Lab-Lift by screwing the rod into the pre-drilled hole.

### 17" Support Rod Kit includes:

- 1 17" (432 mm) Stainless Steel Rod
- 1 Jam Nut
- 1 Flat Washer

Description	Rod Diameter	Item Number
17" Support Rod Kit	0.51" (13 mm)	30400050

#### **Ratchet Tool**

Designed to add extra leverage to your  $12 \times 12^{\circ}$  or  $16 \times 16^{\circ}$  Lab-Lift. This recommended Ratchet Tool easily attaches to the actuating rod to allow easy, accurate adjustments.

Description	Item Number
Ratchet Tool	30400049





### **OHAUS Corporation**

Headquartered in Parsippany, NJ, OHAUS Corporation manufactures an extensive line of high-precision electronic and mechanical balances and scales that meet the weighing needs of virtually every industry. We are a global leader in the laboratory, industrial, and education channels as well as a host of specialty markets, including the food preparation, pharmacy and jewelry industries. An ISO 9001:2008 manufacturer, OHAUS products are precise, reliable and affordable, and are backed by industry-leading customer support. For additional information, contact OHAUS Corporation at 1.973.377.9000, or visit www.ohaus.com.

#### **OHAUS CORPORATION**

\*7 Campus Drive Suite 310 Parsippany, NJ 07054 USA

Tel: 800.672.7722 973.377.9000 Fax:973.944.7177

#### www.ohaus.com

With offices throughout Europe, Asia, and Latin America

\*ISO 9001:2008 Registered Quality Management System